

PRISM CPS Submission Information

PRISM CPS screens are standard offerings for DMSO-soluble compounds in a combination format. These screens are run using our standard design: 5-day assay, 900+ cell lines, test agents run at 7 pt. dose for combinations, with 3-fold dilutions in triplicate. Collaborators select the top screening dose as a part of the screen submission process. Additional workflow information can be found on the [PRISM website](#)

Test Agent Requirements:

- We require at least 400µL per test agent per slot at 1000X the top screening dose in 100% DMSO.
 - For one combination, you will need to submit at least 800µL of your treatment test agent and at least 400µL of your anchor test agent
- All CPS submissions require at least one combination and associated treatment test agent alone (minimum two slots)

Slot #	Treatment Test Agent (7 pt. Dose, 3x dilution)	Volume	Anchor Test Agent (1 dose)	Volume
1	Drug A	400µL	n/a	n/a
2	Drug A	400µL	Drug B	400µL
3	Drug A	400µL	Drug C	400µL

- *Example:* For a drug combination of Drug A +/- Drugs B & C, you would need to submit:
 - Drug A → 1200µL
 - Drug B → 400µL
 - Drug C → 400µL
- For a standard CPS submission, we will screen:
 1. Treatment compound at dose (7-point dose with 3-fold dilutions) alone (Slot #1)
 2. Treatment compound at dose (7-point dose with 3-fold dilutions) + Anchor compound at one anchor dose (Slot #2)
 3. Anchor compound at one anchor dose alone (no slot - since it just 1 dose)
- The collaborator is responsible for ensuring that the test agent is soluble in 100% DMSO at the stock concentration and for QC'ing test agents ahead of submission.
 - We do not have the bandwidth to do any optimization for, or re-running of, test agents that do not perform well in PRISM. Any test agents submitted will be run 'at-risk'
 - We are unable to solubilize powdered stock of compounds for collaborators.