

PRISM EPS Submission Information

PRISM EPS screens are a 10-day assay offering for DMSO-soluble compounds only in a single-agent format. These screens are run using our EPS design: 10-day assay, 500+ solid tumor cell lines, test agents run at 5 pt. dose, with custom dilutions in triplicate. Test agents are re-dosed on Day 6. Collaborators select the top screening dose and dilution factor as a part of the screen submission process. Additional workflow information can be found on the [PRISM website](#)

Test Agent Requirements:

- We require at least 600 μ L at 1000X the top screening dose in 100% DMSO.
 - i.e. If your top screening dose is 10 μ M, you will submit a 10mM stock.
 - The required volume is dependent on the dilution factor the collaborator selects in the submission form.
 - 3-fold dilutions or higher require 600 μ L.
 - 2 to 3-fold dilutions require 720 μ L.
- The collaborator is responsible for ensuring that the test agent is soluble in 100% DMSO at the stock concentration and for QC'ing test agents ahead of submission.
 - We do not have the bandwidth to do any optimization for, or re-running of, test agents that do not perform well in PRISM. Any test agents submitted will be run 'at-risk'
 - We are unable to solubilize powdered stock of compounds for collaborators.